


HÖRLUCHS®



MADE IN GERMANY



ACCESSORIES



HEAR. PROTECT. ENJOY!

CONTENT

HÖRLUCHS® ACCESSORIES

PAGE

About Hörluchs®	3
Hörluchs® Onlineshop	4
Impression material	5
Impression accessories	6
Earmould accessories	8
Milling tools	11
ICP accessories	14
Hearing Protection accessories	15
Filter technology	16
In-Ear accessories	18
Basic equipment	19

Hörluchs® focuses on special solutions at three locations in Bavaria. Whether earmoulds, hearing protection or in-ears – the Hörluchs® logo guarantees quality, comfort and safety. Made in Germany.



HEAD OFFICE

91217 Hersbruck (Ellenbach)

- Administration
- Customer Service
- Field Sales
- Marketing
- Development
- ca. 800 qm



PRODUCTION

91217 Hersbruck (Weiher)


- Order entry
- Audiology
- Modelling
- Production SE / HP / IEM
- Dispatch and warehouse
- ca. 2.000 qm



PRODUCTION

97080 Würzburg

- Modelling
- Production SE
- Dispatch and warehouse
- ca. 400 qm



SHOP.HOERLUCHS.COM

HÖRLUCHS® ONLINESHOP

With us, you can conveniently order any product online - either on our website shop.hoerluchs.com or through your company software (ASEGO, IPN, Amparex).

YOUR ADVANTAGES:

- **Reliable:** Order data in your own company software.
- **Fast:** No need for clarification, faster processing in our company.
- **Convenient:** Ordering possible around the clock.
- **Shorten delivery times by scanning impressions yourself:** Add your own scans of the impressions to your online order. Save up to two days of production time and postage costs.
- **Simple:** trigger reproductions and complaints digitally.

IMPRESSION MATERIAL



HL SOFT 50 ML
AUTOMIX 1: ART. 352
AUTOMIX 2: ART. 3861

Colour: blue/white
 Final hardness: 37 shore A
 Ear retention time: 3:00 min.
 Features: universally usable

Unit: 8 Pieces



HL SPEED 50 ML
ART. 3706

Colour: wine red/white → red
 Final hardness: 30 shore A
 Ear retention time: 1:30 min.

New! Colour change during curing, extra fast!

Unit: 8 Pieces



HL MIX 1.000 ML
ART. 3708

Colour: pink
 Final hardness: 35 shore A
 Mixing time: 1:30 min.
 Ear retention time: 2:30 min.

2 x 500 ml
 Base + Catalysator
 2 dosing spoons pink/white

Unit: 1 Set



HL FLOW 50 ML
ART. 3707

Colour: yellow/white
 Final hardness: 25 shore A
 Ear retention time: 2:00 min.

New! ideal for children's ears, especially for DIC fitting

Unit: 8 Pieces



MIXING CANNULA
ART. 353

For Automix 1 double cannulas
 8 mm diameter

Unit: 100 Pieces



CANNULA TIP
ART. 929

For Automix 1 cannulas

Unit: 50 Pieces



T-MIXING CANNULA
YELLOW ART. 3709

For Automix 2 double cannulas
 4,2 mm diameter

Unit: 25 Pieces



T-MIXING CANNULA PINK
ART. 3710

For Automix 2 double cannulas
 5,4 mm diameter


Unit: 25 Pieces



DOSING SPOON
ART. 3784

For A-Silicones
 1x pink, 1x white

Unit: 1 Pair



MIXING BOARD
ART. 3778

With circular dosing help

Unit: Piece



MIXING SPATULA
ART. 3796

With wooden handle

Unit: Piece



IMPRESSION GUN
ART. 3831

For impression material HL MIX

Unit: Piece

IMPRESSION ACCESSORIES



IMPRESSION PADS 7 MM
ART. 3792

Green, with ventilation

Unit: 24 Pieces



IMPRESSION PADS 9 MM
ART. 3793

Black, with ventilation

Unit: 24 Pieces



IMPRESSION PADS 11 MM
ART. 3794

Beige, with ventilation

Unit: 24 Pieces



IMPRESSION PADS 13 MM
ART. 3795

Blue, with ventilation

Unit: 24 Pieces



IMPRESSION PADS 8 MM
ART. 558

cylindrical

Unit: 100 Pieces



IMPRESSION PADS 12 MM
ART. 559

cylindrical

Unit: 100 Pieces



IMPRESSION PADS 12 MM
ART. 354

conical

Unit: 100 Pieces



IMPRESSION FOIL
ART. 600

For Impressions

Unit: 10 Pieces



BITEBLOCK
ART. 2670

For Impressions

Unit: 5 Pieces



ORDER BAG
ART. 782

With Hörluuchs® Logo

Unit: 50 Pieces





HEINE OTOSCOPE
ART. 517

Heine mini 3000 black
incl. single-use funnel
(5 x 2,5 mm + 5 x 4,0 mm)


Unit: Piece



SPARE BULB
ART. 518

for Heine otoscope
mini 3000


Unit: Piece



HEINE FUNNEL 2,5 MM
ART. 528

Single-use
2,5 mm diameter

Unit: 10 Pieces



HEINE FUNNEL 4,0 MM
ART. 529

Single-use
4,0 mm diameter

Unit: 10 Pieces



HEINE LIGHT STICK
ART. 519

Heine mini C clip lamp
with acrylic tip

Unit: Piece



SPARE BULB
ART. 520

For Heine Light stick
mini C clip lamp

Unit: Piece



ACRYLIC TIP
ART. 1025

For Heine Light stick

Unit: Piece



WAX KNIFE & WAX
ART. 4112

Simply repair defective
impressions with wax

Unit: Set



BAJONET TWEEZERS
ART. 522

Unit: Piece



SCISSORS
ART. 523

Scissors with rounded tip,
for hair removal inside the ear

Unit: Piece



INJECTOR DS 53
ART. 3711

For Automix 2 double cartridges
like HL SOFT (Art. 3861),
HL SPEED and HL FLOW

Unit: Piece



INJECTOR DMA 50
ART. 1026

For Automix 1 double cartridges
like HL SOFT (Art. 352)

Unit: Piece



IMPRESSION CASE SMALL
ART. 683

For home visits incl. 1x Injector pistol, 20x Impression foil, 100x Mixing
cannulas, 5x Funnel 2,5, 5x Funnel 4,0, 1x Scissors, 1x Bajonet tweezers,
100x conical impression pads, 10x Impression material, 50x Pressure seal bag,
20x Order form GS, 20x Function test, 1x auriFIX cleaning spray,
1x Disinfection spray, 1x Otoscope, 1x Light stick, 1x Kitchen paper,
5x Biteblock

Unit: Piece

EARMOULD ACCESSORIES



SOUND TUBE
ART. 3446

Sound tube **2,0 x 3,1 mm** prebent transparent with threading help

Unit: 100 Pieces



SOUND TUBE
ART. 3682

Sound tube **2,0 x 3,1 mm** prebent skin-coloured with threading help

Unit: 50 Pieces



SOUND TUBE
ART. 3256

Sound tube for children **1,5 x 2,5 mm** prebent transparent with threading help

Unit: 25 Pieces



SOUND TUBE
ART. 3280

Sound tube **2,0 x 3,3 mm** prebent transparent

Unit: 25 Pieces



HORN TUBE
ART. 1242

Horn tube **3,0 mm** thin transparent

Unit: 10 Pieces



HORN TUBE
ART. 3687

Horn tube **3,0 mm** thick transparent

Unit: 10 Pieces



HORN TUBE
ART. 1137

Horn tube **4,0 mm** thin transparent

Unit: 10 Pieces



HORN TUBE
ART. 1136

Horn tube **4,0 mm** thick transparent

Unit: 10 Pieces



SOUND TUBE
ART. 686

2,0 x 3,1 mm transparent

Unit: 10 metres



SOUND TUBE
ART. 1195

Sound tube for children **1,5 x 2,5 mm** transparent

Unit: 10 metres



SOUND TUBE
ART. 227

Sound tube **2,0 x 4,0 mm** transparent

Unit: 10 metres



SOUND TUBE
ART. 687

Sound tube stay dry **2,0 x 3,1 mm** transparent


Unit: 10 metres



SOUND TUBE
ART. 1080

Sound tube **2,0 x 3,1 mm** skin-coloured


Unit: 10 metres



CURVED HOSE CONNECTOR
ART. 1174

Curved hose connector **1,5 x 2,5 mm** nylon

Unit: 25 Pieces



CURVED HOSE CONNECTOR
ART. 3705

Curved hose connector **1,5 x 3,0 mm** nylon

Unit: 25 Pieces



HL-TUBELOCK
ART. 3544

Special for silicone earmoulds without gluing

Unit: 10 Pieces




AURIFIX DISPLAY
ART. 3704

Display
incl. 20 auriFIX cleaning sprays

- ✓ Ideal for counter sales
- ✓ With listing of all features
- ✓ Modern design
- ✓ WxHxD: 18,2 x 11 x 14,5 cm

Unit: 1 Set with 20 pieces



AURIFIX COMPRESSED GAS SPRAY
ART. 4310

Unit: Piece



AURIFIX CLEANSING TISSUES
ART. 4311 30 x individually packed
ART. 4312 100 x in dispenser


Unit: Piece



SUPERGLUE
ART. 2913

Special glue 10 g for 3D SOFT SILICONE, 3D TEC LPH & 3D TITANIUM

Unit: Pieces



PRIMER
ART. 303

Primer for superglue 10 ml for 3D SOFT SILICONE Bottle brush


Unit: Piece



SHELLAC
ART. 3835

Shellac for 3D TEC 5 ml bottle brush

Unit: Piece



DRY AIR BLOWER
ART. 3721

For cleaning sound tubes and canals


Unit: Piece



EJECTOR PIN
ART. 2806

Clic mould mini receiver ejector pin

Unit: Piece



DRAWSTRING
ART. 2519

Drawstring

Unit: Piece



MULTIFIX ADAPTER
ART. 3255

Multifix adapter silicone for UNIRIC TEC/TITANIUM


Unit: 10 Pieces



TUBE ANGLE
ART. 2365

Sound tube angle set for 3D SOFT SILICONE (CAST)


Unit: 2 Pieces



HOT AIR TOOL TYPE 25
ART. 3842

PASSAT Hot air tool type 25 220 watt heating capacity not controllable


Unit: Piece



HOT AIR TOOL TYPE 32
ART. 3843


PASSAT Hot air tool type 32 400 watt heating capacity controllable (80°-280° C)

Unit: Piece



TUBE EXPANDING PLIER
ART. 4159

Unit: Piece



STETHOCLIP
ART. 4128

Unit: Piece

EARMOULDS ACCESSORIES



COLOUR SAMPLE TITANIUM
ART. 4226

With seven high-quality
TITANIUM earmoulds

Unit: Piece



SLEEVE AS
ART. 1266

Sleeve for Audio Service, Signia,
Hansaton 5r/5b
modified for UNIRIC SOFT

Unit: 10 Pieces



SLEEVE WIDEX
ART. 1259

Sleeve for Widex &
Easy Wear 5r/5b


Unit: 10 Pieces



SLEEVE FLEXFIT
ART. 1260

for Phonak, Unitron, Alta 60/85,
GN Resound (Größe S)

Unit: 10 Pieces



HÜLSEN OTICON

ART. 1256
for Oticon Agil/Standard receiver 5r/5b

ART. 1255
for Oticon Minifit (Gummi)-
receiver, 5r/5b

Unit: 10 Pieces



SLEEVE CONNEXX
ART. 3249

Sleeve for Connexx special for
Thermotec 5r/5b

Unit: 10 Pieces



SLEEVE PHONAK
ART. 3353

for Sonova 4.0 receiver

Unit: 10 Pieces



SLEEVE SURE FIT
ART. 2735

Sleeve for Sure Fit LP new

Unit: 10 Pieces



SLEEVE AS
ART. 2368

3D TEC LPH sleeve for Audio
Service, Sivantos, Hansaton 5r/5b

Unit: 10 Pieces



**HÖRLUCHS®
SMART SOCKET 2.0**
ART. 4411

Incl. admission for HF4 Filter

Unit: 10 Pieces



**HÖRLUCHS®
SMART SOCKET 2.0**
ART. 4412

Incl. admission for ProWax Filter

Unit: 10 Pieces



**HÖRLUCHS®
SMART SOCKET TOOL**
ART. 4426

For placing and gluing
SMART SOCKET 2.0


Unit: Pieces



CERUSTOP
ART. 1057

Filter for protection against
Cerumen

Unit: 1x Blister with 8 Pieces



CERUSTOP SOCKET
ART. 3389

Socket for cerustop filter


Unit: 10 Pieces



PROWAX
ART. 3773

Filter for protection against
cerumen

Unit: 1x Blister with 6 Pieces



HF4 CERUMENFILTER
ART. 1328-B/1328-R

HF 4 Filter BLACK Spender à 16
Filter rot/blau inkl. Wechselstift

VE 1x Spender

MILLING TOOLS



TWIST DRILL Ø 0,8 MM
ART. 261

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH

Unit: Piece



TWIST DRILL Ø 1,0 MM
ART. 262

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH

Unit: Piece



TWIST DRILL Ø 1,2 MM
ART. 263

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH

Unit: Piece



TWIST DRILL Ø 1,4 MM
ART. 264

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH

Unit: Piece



TWIST DRILL Ø 1,6 MM
ART. 265

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH

Unit: Piece



TWIST DRILL Ø 2,0 MM
ART. 266

Special twist drill for additional drilling, cylindrical, for 3D SOFT SILICONE and 3D TEC LPH


Unit: Piece



ROSEHEAD DRILL Ø 2,5 MM
ART. 3713

For Sound tube drilling for 3D SOFT SILICONE and 3D TEC LPH


Unit: Piece



ROSEHEAD DRILL Ø 2,9 MM
ART. 3712

For Sound tube drilling for 3D SOFT SILICONE and 3D TEC LPH


Unit: Piece



HM-MILL Ø 4,0 MM
ART. 253

Pointed

Unit: Piece



HM-MILL Ø 5,0 MM
ART. 279

Ball

Unit: Piece



HM-MILL Ø 6,0 MM
ART. 2387

Cylindrical

Unit: Piece



HM-MILL Ø 6,0 MM
ART. 254

Cone

Unit: Piece



HM-MILL Ø 8,0 MM
ART. 2425


Ball

Unit: Piece



MANDRELL Ø 8,0 MM
ART. 585


Unit: Piece



HM-MILL Ø 4,5 MM
ART. 4126

Pointed, **for left-handed**

Unit: Piece



HM-MILL Ø 6,0 MM
ART. 4125

Cone, **for left-handed**

Unit: Piece

MILLING TOOLS



SANDPAPER HOLDER
ART. 3776 CONICAL
ART. 3798 CYLINDRICAL

Unit: Piece



SANDING LINE
ART. 244

Tear resistant, durable, grain size 320 for earmoulds 3D SOFT SILICONE and 3D TEC

Unit: Metre



SPONGE STONE
ART. 3714

For editing impressions

Unit: Piece



SANDING SLEEVE HOLDER
ART. 252

0,5 x 11 mm for 3D SOFT SILICONE

Unit: Piece



SANDING SLEEVE
ART. 201

0,5 x 11 mm grain size 120 for 3D SOFT SILICONE

Unit: 10 Pieces



SANDING SLEEVE
ART. 3724

0,5 x 11 mm grain size 140 for 3D SOFT SILICONE

Unit: 10 Pieces



POLISHER FOR TEC/TITANIUM
ART. 3689

Topdent orange wheel pre-polish incl. polishing paste

Unit: Piece



POLISHER FOR TEC/TITANIUM
ART. 3690

Topdent pink wheel high-gloss incl. polishing paste

Unit: Piece



HM-MILL TITANIUM
ART. 3816

Ø 6,0 mm HM mill for Titanium Cone

Unit: Piece



SPIRAL MANDRELL
ART. 3817

For polisher Titanium

Unit: Piece



TITANIUM POLISHER
ART. 3818

Pre-polish for Titanium


Unit: 5 Pieces



TITANIUM POLISHER
ART. 3819

Fine polish for Titanium

Unit: 5 Pieces



HM-MILL TITANIUM
ART. 3815

Ø 4,0 mm Mill for Titanium Pointed


Unit: Piece



TWIST DRILL
ART. 3857

Ø 1,0 mm HM Twist drill for Titanium


Unit: Piece



TWIST DRILL
ART. 3858

Ø 1,2 mm HM Twist drill for Titanium

Unit: Piece



TWIST DRILL
ART. 3859

Ø 1,4 mm HM Twist drill for Titanium

Unit: Piece

SET



DRILLING SET ART. 3846

Drilling basic set inclusive:

- Twist drill 0,8 mm
- Spiral drill 1,0 mm
- Twist drill 1,2 mm
- Twist drill 1,4 mm
- Twist drill 1,6 mm
- Twist drill 2,0 mm
- Rosehead Drill 2,5 mm
- Rosehead Drill 2,9 mm
- HM-Mill 4,0 mm pointed
- HM-Mill 5,0 mm ball
- HM-Fräser 6,0 mm pointed
- Sandpaper holder
- Sanding sleeve holder

Unit: Set

SET



DRILLING/MILLING TITANIUM ART. 3864

Drill and mill for editing titanium earmoulds included:

- HM-Mill Titanium Ø 6,0 mm
- Twist drill Ø 1,4 mm
- Spiral mandrell
- 5x Titanium polisher (Pre-polish)
- 5x Titanium fine polish

Unit: Set



MILL B200 ART. 4107

Mill, battery-powered

Unit: Piece



MILL SAESHIN OZ-BLACK-100A ART. 4340

Mill with table control + Forte 100 a

Unit: Piece



MILL SAESHIN FORTE 200 ALPHA ART. 4116

Mill with ceramic bearings

Unit: Piece



MILL SAESHIN TRAUS STRONG B135 ART. 4120

Mill

Unit: Piece



TESTING EARMOULD ART. 3787 (A), 3788 (B), 3789 (C)

Easy (A), middle (B), hard (C)

- ✓ Testing earmoulds Acryl
- ✓ Ideal for exam preparation
- ✓ Different difficulty levels
- ✓ For left and right ear

Unit: 2 Pieces (L/R)



DRILL HOLDER ART. 3855

Holder for Drills and mills

Unit: Piece



DISINFECTION CONCENTRATE ART. 4417

For disinfection of medical instruments (1 liter)

Unit: Piece

ICP ACCESSORIES



SOUND TUBE ANGLE
ART. 2356

For SOWEI ICP

Unit: 2 Pieces



SOUND TUBE ADAPTER
ART. 2353

For SOWEI ICP

Unit: 2 Pieces



CLOSURE UNIT
ART. 2352 for SOWEI ICP and
for size S, M & P sound tube
connections and receivers
ART. 4127
for ITE Hearing system ICP TIK®

Unit: 2 Pieces



CLOSURE UNIT
ART. 2366

For sound tube angle
SOWEI ICP


Unit: 2 Pieces



SOUND TUBE ADAPTER
ART. 4520

For THERMOTec® ICP

Unit: 2 Pieces



CLOSURE UNIT
ART. 4497

Closure unit for
additional hole THERMOTec® ICP


Unit: 2 Pieces



CLOSURE PLUG
ART. 4498

Closure plug for sound channel
THERMOTec® ICP;
for sound tube connections and
receivers size S & M

Unit: 2 Pieces



CLOSURE PLUG
ART. 4525

Closure plug for sound channel
THERMOTec® ICP for receiver
size P

Unit: 2 Pieces



PROBE TUBE
ART. 754
Insitu probe tube standard
ART. 2347
Insitu probe tube
for SOWEI ICP incl. metal tip


Unit: 2 Pieces



HEARING ANGLE
ART. 3681

For AS Hörluks® BTE
and BTE P G5

Unit: Piece



MICROPHONE COVER
ART. 4092

For ITE hearing aid ICP TIK®

Unit: Piece



ITE CASE
ART. 4077

For ITE hearing aid ICP TIK®

Unit: Piece



BATTERY HATCH
ART. 4086

For ITE hearing aid ICP TIK®

Unit: Piece



SPARE CASES
ART. 3580-1 /-2 /-3

Spare case for BTE 16 G5
dark grey/grey/silver

Unit: Piece



SPARE CASES
ART. 3581-1/-2/-3

Spare case for BTE 16 P G5
dark grey/grey/silver

Unit: Piece




SPARE CASES
ART. 4406-1/-2/-3

Spare case for RIC 16 G7
dark grey/grey/silver

VE Piece

HEARING PROTECTION ACCESSORIES



CARRYING STRAP SILICONE
ART. 115

Silicone strap, clothes clip, cord lock red/blue

Unit: Piece



CARRYING STRAP
ART. 698

For professional hearing protection incl. cord lock

Unit: Piece



STRAP DETECTION
ART. 364

Detectable strap, clothes clip, cord lock blue/blue

Unit: Piece



CARRYING STRAP TEXTILE
ART. 3266

Textile strap, clothes clip, cord lock red/blue

Unit: Piece



CLOTHES CLIP
ART. 116

Clothes clip plastic

Unit: Piece



CLEANING LOOP
ART. 120

For cleaning earmoulds

Unit: Piece



CLEANING BRUSH
ART. 548

For cleaning earmoulds

Unit: Piece



CARE GEL AUDINELL 5 ML
ART. 310

Gel for easier insertion for silicone earmoulds

Unit: Piece



HP ETUI L
ART. 3049

With belt loop

Unit: Piece



HP ETUI S
ART. 3178

With belt loop

Unit: Piece



SPIRAL TUBE
ART. 514

With collar clip for SOWEI COM

Unit: Piece



AURE-X HEADPHONES
ART. 3782

Headphones for functional check, special for our digital function system

Unit: Piece



ANGLE ADAPTER
ART. 2638

For Sound tube connection SOWEI COM transparent

Unit: 2 Pieces



LOCKING UNIT
ART. 963-TR

For SOWEI COM with grip transparent

Unit: 2 Pieces



HF 4 CHANGING PIN
ART. 1029

Tool for changing HF 4 Filter

Unit: Piece



COLOUR SAMPLE SILICONE
ART. 2808

35, 40, 60 shore

Unit: Piece

CARBON FILTER

HLFs Filter for SOWEI and SOWEI COM




HLFs1
ART. 111

SNR 28 dB
Unit: Piece



HLFs2
ART. 112

SNR 26 dB
Unit: Piece



HLFs3
ART. 153

SNR 23 dB
Unit: Piece




HLFs4
ART. 836

SNR 25 dB
Unit: Piece



HLFs5
ART. 837

SNR 22 dB
Unit: Piece



HLFs6
ART. 002

SNR 15 dB
Unit: Piece

HLFs Filter for HA and HAWELxs



HLFs1
ART. 111

SNR 32 dB
Unit: Piece



HLFs2
ART. 112

SNR 28 dB
Unit: Piece



HLFs3
ART. 153

SNR 24 dB
Unit: Piece



HLFs6
ART. 002

SNR 15 dB
Unit: Piece

HLF Filter for HAWEI and HAWELs



HLF1
ART. 109

SNR 28 dB
Unit: Piece



HLF2
ART. 110

SNR 26 dB
Unit: Piece



HLF3
ART. 152

SNR 23 dB
Unit: Piece



HLF6
ART. 948

SNR 15 dB
Unit: Piece

HLFxs Filter for SOWEI FREE



HLFxs26
ART. 2697

SNR 26 dB
Unit: Piece




HLFxs20
ART. 2698

SNR 20 dB
Unit: Piece



HLFxs10
ART. 2700

SNR 10 dB
Unit: Piece



HLM FILTER



HLM30
ART. 2922

SNR 30 dB

Unit: Piece



HLM25
ART. 2921

SNR 25 dB

Unit: Piece



HLM20
ART. 2920

SNR 20 dB


Unit: Piece



HLM15
ART. 2919

SNR 15 dB

Unit: Piece



HLM10
ART. 3627

SNR 10 dB

Unit: Piece



HLM SLEEVE
ART. 3661

Für HA Unlimited

Unit: Piece

Hörluchs® has developed a filter system which in combination with the appropriate earmould offers the customer a high degree of choice. In combination with **suitable Hörluchs® hearing protection earmoulds**, the HLM exchangeable filters have a **type examination certificate**. This means that this hearing protection with its filter technology can be used in both the private and commercial sectors.

With the five linear HLM spin-on filters, the ambient noise can be controlled as required. The robust HLM filters from Hörluchs® are easy to replace and durable. In contrast to comparable filters from other manufacturers, they can be easily cleaned with an earpiece in an ultrasonic bath.

CHANGING FILTER



CHANGING TOOL
ART. 612

Filter setting and removal tool for SOWEI

Unit: Piece



CHANGING TOOL
ART. 3052

Filter setting and removal tool for HLM filter

Unit: Piece



CHANGING TOOL
ART. 138

Filter setting and removal tool for HAWEI/s

Unit: Piece



FLAT ROUND NOSE PLIERS
ART. 916

Chromed

Unit: Piece



CHANGING TOOL
ART. 187

Filter setting and removal tool metal for HCP earmoulds

Unit: Piece

IN-EAR ACCESSORIES



CLEANING CLOTH
ART. 2905

Cleaning cloth black

Unit: Piece



SILICONE DOMES
ART. 2641 / 2642 / 2643

Silicone domes for Hörluchs® In-Ears (S/M/L)

Unit: 10 Pieces (of one size)



MEMORY FOAM DOMES
ART. 2644 / 2645 / 2647

Memory foam domes for Hörluchs® In-Ears (S/M/L)

Unit: 10 Pieces (of one size)



PREMIUM CABLE
ART. 2707

Premium cable, black incl. bending element, 1,6 m, changeable, 2 Pin plug

Unit: Piece



PREMIUM CABLE
ART. 3116

Premium cable, black, without bending element, 1,6 m, changeable, 2-Pin plug

Unit: Piece



PREMIUM CABLE
ART. 2706

Premium cable, transparent incl. bending element, 1,6 m, changeable, 2 pin plug

Unit: Piece



EASYUP CABLE
ART. 3193

Cable for EasyUp, black, 1,5 m, changeable, mmcx connection

Unit: Piece



HIGH END CABLE
ART. 3050

High End, black, 1,6 m, changeable, T2 plug

Unit: Piece



LYNX CASE
ART. 2374

Case for In-Ears and HA Active

Unit: Piece



IN-EAR TEST CASE
ART. 4219

Four HL4 testing In-Ears in a handy hardcase

Unit: Piece



IN-EAR AUDIO TOWER
ART. 3639

Audio tower incl. Equipment

- ✓ Professional presentation
- ✓ Comparison of 2 to 6 ways
- ✓ incl. domes & catalogues
- ✓ WxHxD: 49 x 123 x 40 cm

Tablet and songs are not included!

Unit: Piece



BASIC EQUIPMENT SET

THE PERFECT STARTER KIT FOR YOUR BRANCH

Additional **10% discount** on all items in this set*.

ART.	DESCRIPTION	QUANTITY		
1026	INJECTOR DMA 50	1 Piece	✓	
352	HL SOFT 50 ML	8 Pieces	✓	
353	MIXING CANNULA	100 Pieces	✓	
3711	INJECTOR DS 53	1 Piece	✓	
3706	HL SPEED 50 ML	8 Pieces	✓	
3707	HL FLOW 50 ML	8 Pieces	✓	
3709	T-MIXING CANNULA YELLOW	25 Pieces	✓	
3710	T-MIXING CANNULA PINK	25 Pieces	✓	
3708	HL MIX 1000 ML	1 Set	✓	
3778	MIXING BOARD	1 Piece	✓	
3784	DOSING SPOON	2 Pieces	✓	
3796	MIXING SPATULA	1 Piece	✓	
558	IMPRESSION PADS 8 MM	100 Pieces		✓
559	IMPRESSION PADS 12 MM	100 Pieces		✓
354	IMPRESSION PADS 12 MM	100 Pieces		✓
600	IMPRESSION FOIL	10 Pieces		✓
2670	BITEBLOCK	5 Pieces		✓
782	ORDER BAG	50 Pieces		✓
517	HEINE OTOSCOPE	1 Piece		✓
528	HEINE FUNNEL 2,5 MM	10 Pieces		✓
529	HEINE FUNNEL 4,0 MM	10 Pieces		✓
519	HEINE LIGHT STICK	1 Piece		✓
522	BAJONET TWEEZERS	1 Piece		✓
523	SCISSORS	1 Piece		✓
3842	HOT AIR TOOL TYPE 25	1 Piece	✓	
3446	SOUND TUBE	100 Pieces		✓
1174	CURVED HOSE CONNECTOR	25 Pieces		✓
2913	SUPERGLUE	1 Piece		✓
3835	SHELLAC	1 Stück		✓

OPTIONAL
FIX

ART. 4021

*Illustration without optional items



ART.	BEZEICHNUNG	QUANTITY		
3721	DRY AIR BLOWER	1 Piece		✓
261	TWIST DRILL 0,8 MM	1 Piece		✓
262	TWIST DRILL 1,0 MM	1 Piece		✓
263	TWIST DRILL 1,2 MM	1 Piece		✓
264	TWIST DRILL 1,4 MM	1 Piece		✓
265	TWIST DRILL 1,6 MM	1 Piece		✓
266	TWIST DRILL 2,0 MM	1 Piece		✓
3713	ROSEHEAD DRILL 2,5 MM	1 Piece		✓
3712	ROSEHEAD DRILL 2,9 MM	1 Piece		✓
253	HM-MILL 4,0 MM	1 Piece		✓
279	HM-MILL 5,0 MM	1 Piece		✓
2387	HM-MILL 6,0 MM	1 Piece		✓
254	HM-MILL 6,0 MM	1 Piece		✓
2425	HM-MILL 8,0 MM	1 Piece		✓
585	MANDRELL 8,0 MM	1 Piece		✓
252	SANDING SLEEVE HOLDER	1 Piece		✓
201	SANDING SLEEVE	10 Pieces		✓
3798	SANDPAPER HOLDER	1 Piece		✓
244	SANDING LINE	5 Metres		✓
3714	SPONGE STONE	1 Piece		✓

OPTIONAL
FIX

*A 10% discount applies to all items in this set. Optional items can be added individually.

HEAD OFFICE & INVOICE ADDRESS

Hörluchs Hearing GmbH & Co. KG
Bergseestraße 10 · 91217 Hersbruck · Germany
Tel. +49 (0)9151 90886-0
Fax +49 (0)9151 90886-27



PRODUCTION & ORDER ENTRY

Friedrich-Köchert-Str. 4 · 91217 Hersbruck · Germany

PRODUCTION

Würzburg · Germany

CUSTOMER SERVICE info@hoerluchs.com

ICP / AUDIOLOGIE audiologie@hoerluchs.com

MARKETING marketing@hoerluchs.com

Our GTC apply in the currently valid version. Errors, misprints or other mistakes are reserved and Hörluchs® assumes no liability. V.i.S.d.P. Thomas Meyer, Managing Director & Owner. Item no.: 4522 • Version: 20240119a



HÖRLUCHS®